


## Hand Held Refractometer MASTER-53T


 Shop

 Quotation Request

 Send an inquiry

MASTER-53T are capable of measuring not only various kinds of juice, coffee, and other beverages but also processed food such as sauce, ketchup, low calorie jam, etc. These are an Automatic Temperature Compensation (ATC) model. The MASTER-T / PT series is a good alternative choice when the [water resistant feature](#) is not a priority. The MASTER-53T is made of metal, and the MASTER-53PT is plastic. The specifications are exactly the same. For salty and acidic samples, choose the plastic body, MASTER-53PT.

#### Specifications

Model	MASTER-53T
Cat.No.	2352
Range	Brix : 0.0 to 53.0%
Minimum Scale	0.2%
Accuracy	±0.2%
Repeatability	±0.1%

#### Options

- SUPER EXTRACTOR : RE-29401
- 10% Sucrose (±0.03%) : RE-110010
- 20% Sucrose (±0.03%) : RE-110020
- 30% Sucrose (±0.03%) : RE-110030
- 40% Sucrose (±0.04%) : RE-110040
- 50% Sucrose (±0.05%) : RE-110050
- Daylight plate for small volume samplings : RE-2311-67M
- Daylight plate : RE-2315-60M
- \* Calibration Certificate : Contact an ATAGO representative for details.